



EyeReader

Webcam-based eye tracking technology

By Tess den Uyl and Amogh Gudi

Research Consultants, Vicarious Perception Technologies BV



A solution by Noldus Information Technology

WEBCAM-BASED EYE TRACKING TECHNOLOGY

You can use EyeReader in a flexible and mobile lab environment or within the FaceReader Online platform for online testing.

Next to facial expressions analysis, we now also developed eye tracking analysis, performed with a regular webcam. Our eyes are our windows to the world. Our gaze reveals what we are looking at and this often indicates what we find important. Eye gaze can be stimulus driven and reveal what implicitly draws attention, but is also driven by explicit motivation.

Therefore, eye tracking is used in a wide domain of fields, such as psychology, marketing, and usability research. Usually, eye tracking technology consists of an infrared camera, attached to the computer (with different price ranges and accuracy, depending on the requirements). However, using a simple RGB camera enables eye tracking to be done from every computer/laptop with a camera, thus increasing flexibility greatly.

This white paper discusses the methodology and accuracy behind the technology. It also gives a validation example for a specific application and indicates how to start using it in your own project.




HOW IT WORKS

The EyeReader algorithm estimates the gaze direction and relates that to an image on the screen. The first step is face analysis to locate key-points in the face (such as eyes, nose, etc.) and estimate the head pose.

The eye image and head pose information are then used by a pre-trained neural network (a convolutional neural network or CNN) to produce a 3D gaze vector (a line segment in three directions, x,y,z) representing the direction of the person's eye gaze. The neural network is trained with many labelled data-sets of screen locations and video recordings to learn the relation between the image of the eyes and the gaze vector.

To go from a 3D model of the gaze to a location on a screen one needs to know exactly where the screen is. With the help of a calibration task in which the participant follows dots on a screen, the 3D vector can be mathematically decoded to 2D gaze x,y points on the screen.

Illustration of processing steps of EyeReader




ACCURACY OF EYEREADER

We have tested on our own validation dataset, which includes data from 17 participants (12 males/5 females, 9 aged 18-29/8 aged 30-59) with EyeReader and a Tobii Nano eye tracker with natural lighting conditions. Participants performed a calibration task and a few psychological tasks. They were situated with a distance of approximately 75 cm in front of a 1920x1080 23-inch screen and also performed a calibration task on a 15-inch laptop (without eye tracker).

The results on the 32-dot calibration task (different dots used for calibration and testing) on a desktop computer show that the system predicts gaze points on screen with an average deviation of 2.4 ± 0.6 cm (or an average of 113.6 pixels, and average of 5.2% deviation of the screen).

Also, the results demonstrate that a large amount of sections of the screen can be discerned with high accuracy (it can estimate a gaze location in $1/16^{\text{th}}$ of the screen with 96% accuracy), albeit not in the range of hardware eye tracking systems (having over 99% accuracy for $1/16^{\text{th}}$ of the screen). If one wants to study a difference in fixation between two points on the screen, it is ideal if they are in different $1/16^{\text{th}}$ sections of the screen.

Example of X and Y gaze point estimates of 1 participant for EyeReader and the Tobii Nano eye tracker.




The results are comparable with competitor systems on the market, although it should be noted that direct comparison is not possible due to different screen sizes (making noting errors in pixels problematic) and different testing settings.

LAPTOP VERSUS DESKTOP

Another outcome of this test was that the system is more accurate on the laptop (Mean Absolute Error: 2.3 ± 0.5 cm) in comparison with the desktop (MAE: 2.5 ± 0.7 cm, also see figure below). This is likely due to the smaller screen and the fixed ideal alignment of the in-built webcam of the laptop. The results are similar for brown (MAE: 2.5 cm) in comparison to blue eyes (MAE: 2.4 cm), and for those not wearing glasses (MAE: 2.4 cm) in comparison to those that are wearing glasses (MAE: 2.7 cm). However, this is based on only 3 people who are wearing glasses (with thin lenses); when there is reflection of the surrounding light this will likely cause issues.

Image of gaze point estimation accuracy for the two devices with different tolerance levels (deviation classified as accurate).



HOW CAN EYEREADER BE USED


An important but simple finding within eye tracking research is that people look longer at what they find interesting. We let the participants perform two simple choice tasks, where they viewed a few images of items (e.g. cookies, chocolate, etc.) and were allowed to choose one afterwards.

A comparison of the EyeReader gaze estimations with the eye tracker estimations shows a high correlation between the X and Y coordinate estimation (X: Mean $r = 0.87$, Y: Mean $r = 0.77$).

Gaze estimations can be classified into being a saccade (fast movements with less visual information processing) or fixation (slower eye focus relevant for information processing) based on the speed of the gaze movements. For each category that was presented total duration of fixation time was calculated.


VALIDATION OF THIS APPLICATION

Both systems resulted in the same interpretation of the results: 79% of participants looked longest at the category that they preferred (choice as a reward for participation), with EyeReader giving an average of 9.3 sec on



The image shows a screenshot of a web page titled "Behavioral Research Blog" by Noldus Information Technology. The page features a navigation menu at the top with links for "NOLDUS.COM", "PRODUCTS", "APPLICATIONS", "CUSTOMER STORIES", and "CONTACT". A search bar and a subscription form are located in the top right corner. The main content area displays a blog post titled "Creating a custom expression for Engagement: A validation study with FaceReader". The post is attributed to Deyenne Sarkol and dated Monday, January 2021. The text below the title reads: "By Deyenne Sarkol-Teulings, graduate student at Applied Cognitive Psychology at the University of Utrecht, The Netherlands". The entire page is overlaid with a network of blue dots and lines, representing eye-tracking data points and movement paths across the content.

Gaze results can help track what captures the participant's attention.



Correlation within different categories in two tasks between fixation duration of Tobii Nano eye tracker and EyeReader output.

chosen category and 2.0 sec on not chosen categories (versus 11.1 sec and 2.6 sec for the eye tracker).

The slightly lower estimates in EyeReader are likely due to a slightly higher percentage of failed eye localizations. In addition, the total fixation duration for each category indicated strong correlation between the two systems (Mean $r = 0.91$, see figure above).

HOW TO IMPLEMENT EYEREADER

EyeReader is very suitable for researching realistic behaviors where you expect clear discernible effects.

The results show that eye tracking via the webcam can be a useful and easily applicable research tool. Under good conditions varying individual and screen characteristics do not cause large differences in performance. Generally speaking, there is often a trade-off between ecological validity and measurement precision.

If high precision is required, a lab setting is more appropriate where participants can also put their head in a headrest and have long calibration sessions. EyeReader is very suitable for researching realistic behaviors where you expect clear discernible effects.

You can use EyeReader in a flexible and mobile lab environment or within the FaceReader Online platform for online testing. If you are already a FaceReader Online user, make sure you are using FaceReader 9. When you open a project, click on the tab 'details' and make sure version 9 is selected as analysis engine for the participants' recordings. You can then turn gaze tracking on in your experiment design.

Go to www.facereader-online.com for more information and sample studies. To start using FaceReader Online or for more information on how to integrate EyeReader in a lab environment or other tool, please contact [Noldus](#).



INTERNATIONAL HEADQUARTERS
Noldus Information Technology bv
Wageningen, The Netherlands
Phone: +31-317-473300
Fax: +31-317-424496
E-mail: info@noldus.nl

NORTH AMERICAN HEADQUARTERS
Noldus Information Technology Inc.
Leesburg, VA, USA
Phone: +1-703-771-0440
Toll-free: 1-800-355-9541
Fax: +1-703-771-0441
E-mail: info@noldus.com

REPRESENTATION
We are also represented by a
worldwide network of distributors
and regional offices. Visit our website
for contact information.

WWW.NOLDUS.COM